

Takahe Analytical Instruments B-104, Punit Industrial Premises Plot 11 & 11 A, TTC Industrial Area, MIDC Turbhe Navi Mumbai 400 703 INDIA Tel: +91(22) 2087 4443 | +91(22) 2087 9852

Email: <u>takahe.instruments@gmail.com</u> info@takaheinstruments.com

Website: <u>www.takaheinstruments.com</u>

Speedovap – High Volume -Combo Nitrogen Evaporator & Concentrator

About Us: Takahe Analytical Instruments is a manufacturer & supplier of bio- analytical & quality testing instruments such as Nitrogen Evaporators and Concentration System, Solid Phase Extraction(SPE)-Positive Pressure Manifold besides supplying HPLC spares and Second hand Refurbished & Reconditioned HPLCs.

Nitrogen Evaporators: Nitrogen evaporators are used to prepare samples for analysis. Used in several industries, such as food, biotechnology, analytical chemistry, petroleum and environmental sciences, Nitrogen Evaporators prepare sample by evaporation using combination of heat and nitrogen to concentrate samples for further testing on systems like HPLC, GC or Mass Spectrometry.

Takahe Analytical Instruments an Indian manufacturer of Nitrogen Evaporators & Concentration systems offer a range of models of the Nitrogen Evaporators. Our flagship product has been satisfactorily accepted by end-users. To serve all the industry needs we have also in our offering, a High Volume Nitrogen Evaporator often used for Food, Soil and pesticide testing and other such needs.

With over 600 installations across India and overseas, Takahe Analytical Instruments have established their product with CROs, agriculture & environment industry, food industry, testing labs and educational institutions. While adding instruments to the lab, users after having used the product choose this product again as their choice of evaporators.

Specifications & Technical Details:



Standard Model: Details on right. Customisation for available after feasibility study for users having a requirement of Low Volume and High volume evaporation. (COMBO)

Model name: Speedovap-HV & Combo

Max no. Of samples*

Upto 8 /15/ 50

Sample Volume

2 ml to 250 ml

Number of Stations 2/3 / 5

Max no of samples /station Upto 3/5/10 *

Evaporation Medium Nitrogen

Temperature Setting Ambient to 90°C

Pressure Adjustable

Time setting 1 min to 9.5 hrs or indefinite

*Number of samples dependant on glassware shape & volume



Takahe Analytical Instruments
B-104, Punit Industrial Premises
Plot 11 & 11 A, TTC Industrial Area, MIDC
Turbhe Navi Mumbai 400 703 INDIA
Tel: +91(22) 2087 4443 | +91(22) 2087 9852
Email: takahe.instruments@gmail.com

info@takaheinstruments.com

Website: www.takaheinstruments.com

Speedovap - High Volume Nitrogen Evaporator & Concentrator

Other details:

Controls Microprocessor

Temp & Time display LED for better visibility

Supplied with Installation Kit:

Drain water pump Exhaust hose Nozzle plugs Nozzle cleaner Power cord

Manual, IQ/OQ/PQ docs with traceability certificates

Rating 230 volts

Features:

- Aluminum body for rust control
- In-built alarms with customizable setting for over temperature alert
- Sensor for Door Open triggering alert
- Inter-changeable test tube racks for multiple size test tubes testing
- Safety Features include Audio and display Alarms; Auto Stop of process on Door Open etc.

SPECIAL FEATURE:

- Customisation of racks for different kind of glassware
- End Point test tubes available on specific request for timed end point method

Certifications:

International Certification conforming to

- 2004/22/EC Measuring Instruments
 Directive
- 2014/35/EU Low Voltage Directive
- 2009/105/EC Simple Pressure Vessels

ISO Certification to Takahe Analytical Instruments

 ISO 9001:2015 Quality Management System for following scope of activities: Manufacturing of Analytical Instruments and Related Software Development.



